

Solution to form ideal & aesthetic emergence profile for each tooth position

Anatomic Healing Abutment™

- Facilitates aesthetic final prosthesis, like natural teeth
- Shapes fit perfectly with extraction socket
- Allows immediate loading, reducing chairtime




Optimal gingiva healing

- Various designs to match tooth shape of each extraction socket
- Save time, save money
- Ready-made products, no need for resin or milling



Aesthetic final prosthesis, like natural teeth

Supports same gingiva form as with natural teeth to expedite aesthetic final prosthesis
Satisfactory results for both clinicians and patients



Optimal collaboration with digital abutment

Preservation of natural gingiva form using Anatomic Healing Abutment™ facilitates completion of final patient-customized abutment and prosthesis

Clinical Cases

Courtesy of Dr. Sam Omar

Case 1.



Case 2.




► Anatomic Healing Abutment™

All in one package



Incisor



Incisor				
MD (mm)	4.0	4.5	6.0	7.0
LL (mm)	5.0	4.5	5.0	6.0
Position No.	I1	I2	I3	I4

Canine



Canine	
MD (mm)	5.0
LB (mm)	5.5
Position No.	C1

Pre-Molar




Pre-Molar		
MD (mm)	4.5	5.0
LB (mm)	6.0	7.0
Position No.	P1	P2

Molar



Molar								
MD (mm)	6.0	6.0	6.0	7.0	7.0	7.0	8.0	8.0
LB (mm)	7.0	8.0	9.0	8.0	9.0	10.0	9.0	10.0
Position No.	M1	M2	M3	M4	M5	M6	M7	M8

Special



Special											
MD (mm)	4.5	5.0	5.0	6.0	6.0	6.0	7.0	7.0	7.0	8.0	8.0
LB (mm)	6.0	6.5	7.0	7.0	8.0	9.0	8.0	9.0	10.0	9.0	10.0
Position No.	S1	S2	S3	S4	S5	S6	S7	S8	S9	S10	S11

Order code

System	Height(mm)	Type	Ref.c
ANYRIDGE® by MEGAGEN	4	Hex	C-AHAR04HP
		Non-Hex	C-AHAR04NP
	5	Hex	C-AHAR05HP
		Non-Hex	C-AHAR05NP
	7	Hex	C-AHAR07HP
		Non-Hex	C-AHAR07NP
BLUEDIAMOND IMPLANT	4	Non-Octa	C-AHBD04NP
		Octa	C-AHBD04OP
	5	Non-Octa	C-AHBD05NP
		Octa	C-AHBD05OP
	7	Non-Octa	C-AHBD07NP
		Octa	C-AHBD07OP
ANYONE® by MEGAGEN	4	Hex	C-AHAO04HP
		Non-Hex	C-AHAO04NP
	5	Hex	C-AHAO05HP
		Non-Hex	C-AHAO05NP
	7	Hex	C-AHAO07HP
		Non-Hex	C-AHAO07NP