


Metal-free Prosthetic Solution



AXA Abutment

- Full mouth prosthetics with no extra components
- Advanced Prosthetics Solution makes the impression easier by reducing the chair time
- Choose the best fit from the diverse prosthetic options
- Support various restorative techniques, including Multi case and Cementation




Metal-free

A diverse selection of Angled abutments reduces the chair time by enabling direct connection to crowns and enhances crown fabrication convenience



*FEA for Stress distribution





AXA Abutment

	Multi-unit Abutment	AXA Abutment
Post Height	2.2mm	4mm
Complete Screw	1.5 Thread	3.5 Thread
Screw Strength Enhancement	M1.4	M1.6


Optimal design for crown stress distribution (prevent chipping)

Optimal design for fracture prevention enables to check stress distribution from FEA



Precise and easy fabrication of prosthetics

Compared with the existing Scan Abutment, the prosthetic configuration rate has quadrupled in Multi Case



Truly Cement less

Directly connect the crown and the screw without intraoral and out-oral cementation



Click on connection makes everything simple

Achieve easier attachment and detachment of Temporary Crown with Click on




Abutment Screw per Cuff types

Enhanced accessibility for the Hand Driver


Minimized unnecessary screw hole: Offer a Screw that is dedicated to each cuff

ARI®
by MEGA^{GEN}





Angle	Cuff	
	Max Cuff	Middle Cuff
0°	2	
	3	
	4	
	5	
12°	3	2.5
	4	3.5
	5	4.5
20°	3	2.5
	4	3.5
	5	4.5
30°	3	2
	4	3
	5	4


BLUEDIAMOND IMPLANT



Angle	Cuff	
	Max Cuff	Middle Cuff
0°	2	
	3	
	4	
	5	
12°	3	2.5
	4	3.5
	5	4.5
20°	3	2.5
	4	3.5
	5	4.5
30°	3	2
	4	3
	5	4



Clinical Case 1. Mandibular Full Case



Clinical Case 2. Maxilla Full Case



Clinical Case 3. ARi Bridge Case

